



## TRAVELER

### Value and practicality for basic weighing

The OHAUS Traveler™ Series was designed with a unique draftshield that offers weigh-thru capability while being used and 'stackability' when being stored! With the right combination of performance and features, the Traveler™ is great value for your basic weighing applications.












Construction    ABS plastic with stainless steel pan

Design Features    Unique Weigh-thru draftshield designed with an easy to remove centre piece that provides protection from rough environments while maintaining weighing speed, designed to stack, shipping lock, calibration lock, fixed integral weigh below hook, calibration weight supplied with 150g and 300g models

Display    LCD

Operation    Mains (included), optional battery with programmable Auto Shut-Off

Communication    Optional RS232 or USB interface with cable

Photo	Model Code	Item Nr.	Capacity	Readability	Platform		Test Weight
	TA152	72202482	150g	0.01g	120mm	X	150g, OIML F2 or ASTM Class 4
	TA302	72202483	300g	0.01g	120mm	X	300g, OIML F2 or ASTM Class 4
	TA301	72202484	300g	0.1g	120mm	X	300g, OIML F2 or ASTM Class 4
	TA501	72202485	500g	0.1g	120mm	X	500g, OIML F2 or ASTM Class 4
	TA1501	72202486	1500g	0.1g	120x135mm	X	1500g, OIML F2 or ASTM Class 4
	TA502	83997626	500g	0.01g	120mm	X	500g, OIML F2 or ASTM Class 4
	TA3001	72202487	3000g	0.1g	120x135mm	X	3000g, OIML F2 or ASTM Class 4
	TA5000	72202489	5000g	1g	120x135mm	X	5000g, OIML F2 or ASTM Class 4